

Bⁱ screen

Protective face shields



Protection shield
against splashes
and drops

- Design
- Safety
- Comfort

Bⁱ

Bⁱzeta

thermoplastic extrusions

Bⁱ screen

Protective face shields

ADVANTAGES

- Face shield can be reused indefinitely and easily cleaned and sanitized.

- They are comfortable to wear, reduce the potential for inoculation by preventing the wearer from touching their face.


- No need to be removed to communicate.

- All ways remember to maintain social distancing.

- Allows visibility of facial expression and lip movements for speech perception.

- Colleagues and customers can easily recognise you.

- The product is recyclable and unlike daily disposable masks this will give length of service and more cost effective whilst being kinder to the environment.



Design by **MOREDESIGN**

MADE IN ITALY

bizeta.it

Bizeta s.r.l.
Via delle industrie, 11
30030 Salzano (VE)
info@biscreen.it

This protective device has been
designed in accordance with the
EN:166 2004 specifications

Bⁱzeta

thermoplastic extrusions

Bⁱ screen

Protective face shields

Headband in thermoplastic material

forehead band is
suitable to skin contact,
water and sweat resistant

Adjustment
elastic band

B-SCREEN PRO

Polycarbonate screen, excellent optical properties, great thermal and surface resistance. Headband in material suitable for medical use

BASE VERSION

Polyester screen, resistant to chemicals, any optical distortion. Suitable to food contact

STANDARD COLORS



BLUE



CHARCOAL GREY



WHITE

The colors are purely indicative.

CUSTOMIZATION



BRANDING: corporate brand or image application

SPECIAL COLOURS: on request

bizeta.it

Bizeta s.r.l.
Via delle industrie, 11
30030 Salzano (VE)
info@biscreen.it

This protective device has been
designed in according to the
EN:166 2004 specifications

Bⁱzeta

thermoplastic extrusions

Bⁱscreen

Protective face shields



Instructions

FOREHEAD APPLICATION


1_The forehead section is shaped and equipped with adhesive. Remove the adhesive protection and apply the forehead into the space delimited by two lubber lines in relief on the surface.

2_Fold the final edge of the elastic band to get a coupling on two slots. Proceed to do the same with the other side.

Easy band



Easy Click



bizeta.it

Bizeta s.r.l.
Via delle industrie, 11
30030 Salzano (VE)
info@biscreen.it

This protective device has been
designed in according to the
EN:166 2004 specifications

Bⁱzeta

thermoplastic extrusions